



**KANTECARE**  
we care, wherever you are



**Luer Lock Safety Syringe**

**Item No.:** KTC28-LLSS

### **Product Details**

Kantecare Luer Lock Safety Syringe was engineered from the ground-up over a number of years to be the best safety syringe available. The Luer Lock Safety Syringe provides assured protection and performance advantages for healthcare professionals. The Luer Lock hub means the syringe is compatible with all kinds of standard hypodermic needle.

### **Designed for Safety**

- Engineered, integral safety design
- No add-on pieces
- Hands and fingers stay behind the needle at all time
- Needle withdraws into the barrel
- Permanently disabled

### **Easy to Use**

- Looks and functions like a traditional syringe
- Minimal training and supervision required
- No change in clinical technique
- No obstructing parts
- Optimal low-angle injection
- Audible and visible confirmation of safety features
- Bold, easy-to-read gradations

### **Economical**

- Accurate dose with minimal dead space
- Reduce hazardous waste volume
- Potential to reduce waste disposal cost

**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: info@kantecare.com MSN: kantecare@hotmail.com

SKYPE: kantecare



[www.kantecare.com](http://www.kantecare.com)



# **KANTECARE**

**we care, wherever you are**

- Inventory standardization with full product range
- Ready compliance with industry regulations
- Reduce costs associated with needlesticks



**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: [info@kantecare.com](mailto:info@kantecare.com) MSN: [kantecare@hotmail.com](mailto:kantecare@hotmail.com)

SKYPE: kantecare

**[www.kantecare.com](http://www.kantecare.com)**



**KANTECARE**  
we care, wherever you are



**Without Needle**

**Item No.:** KTC28-LLSS01

#### **Product Details**

Kantecare Luer Lock safety syringe without needle is part of our safety injection device solutions. This application offers the added ability to use all kinds of standard needle to draw up medication or injection. The Luer Lock nozzle also means the syringe is compatible with needle-free connector system. Our Luer Lock safety syringe without needle is available in sizes from 2ml to 20ml.

**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: [info@kantecare.com](mailto:info@kantecare.com) MSN: [kantecare@hotmail.com](mailto:kantecare@hotmail.com)

SKYPE: kantecare



[www.kantecare.com](http://www.kantecare.com)



**KANTECARE**  
we care, wherever you are



**With Detachable Needle**

**Item No.:** KTC28-LLSS02

#### **Product Details**

Kantecare Luer Lock safety syringe with detachable needle is part of our safety injection device solutions. This application offers the added ability to exchange needles or use our cannula (the blunt syringe filling needle and Plastic Vial Access Plastic Cannula) which is used to replace the standard needle to draw up medication from Vials before injection. Our Luer Lock safety syringe with detachable needle is available in sizes from 2ml to 20ml with a wide variety of needle sizes.

**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: [info@kantecare.com](mailto:info@kantecare.com) MSN: [kantecare@hotmail.com](mailto:kantecare@hotmail.com)

SKYPE: kantecare



**[www.kantecare.com](http://www.kantecare.com)**



**KANTECARE**  
we care, wherever you are



### Blunt Fill Needle

Item No.: KTC28-LLSS03

#### Product Details

Kantecare blunt fill needle is used to prepare medication prior to injection as an effective way to reduce needlestick injuries.

Suitable for use with our Luer Lock safety syringes, as well as all standard Luer tip syringes.

#### Features

Effectively reducing needlestick injuries related to vial penetration.

- A lower cost alternative than sharps safety engineered devices.
- Compatible with all standards Luer fit syringes.
- Easy to identify by red color.

#### Caution!

Our blunt fill needle is designed for medication preparation only and cannot be used for injection.

**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: info@kantecare.com MSN: kantecare@hotmail.com

SKYPE: kantecare



[www.kantecare.com](http://www.kantecare.com)



**KANTECARE**  
we care, wherever you are



**Fixed Needle Safety Syringe**

**Item No.:** KTC28-FNSS

### Product Details

Kantecare fixed needle safety syringe has been specifically designed to aid in the prevention of syringe reuse and needlestick injury. New enhanced feature make it the right device for allergy injection, insulin injection, TB injection and both large scale and routine immunization program.

Size available: 0.1ml, 0.5ml, 1ml, 30U, 50U, 100U

### Features

- Minimal Dead space.
- High Quality, permanently attached needle.
- Patented needle retraction mechanism allow user to withdraw the needle back into the barrel.
- Patented Locking ring design allow user to lock the needle into the barrel after needle retracted.
- Patented No Gap\* design.
- Easier to remove air bubble.
- Less wastage compare with standard and other safety syringe design.
- Clear Barrel and Bold scale markings.
- Individual blister packing.
- Low Injection/pull force.
- Minimum training required.
- Affordable price.



**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: info@kantecare.com MSN: kantecare@hotmail.com

SKYPE: kantecare

[www.kantecare.com](http://www.kantecare.com)



# KANTECARE

we care, wherever you are

## Benefits

- Greater Patient comfort
- Minimal Vaccine waste
- High effective in preventing syringe reuse and needle stick injuries
- Take less space in Sharp container; significantly reduce the cost for dispose syringes and needles
- Ease to see and measure syringe content
- Accurate dose measurement
- Environmentally friendly
- Maintains content sterility when package unopened and not damaged
- Improved to access to safety technology

## Minimal dead space

Using a low waste space syringe reduces medication waste and medication costs and allows for accurate dosing. Kantecare fixed needle Safety Syringe has a specially designed needle hub and engagement mechanism. This unique design reduces the amount of waste space volume typically found in devices with a Luer design and other designs of safety syringes. The waste space volume for Kantecare fixed needle safety syringe is 0.008ml which easily meets the standard requirements of a conventional insulin syringe, and 90% less than 1ml syringe with a detachable needle (.07 ml).

## No gap design

Kantecare fixed needle safety syringe offers the same no gap design as a conventional fixed needle syringe, allowing the clinician to draw up medication from the zero graduation mark and easy to remove the air bubbles from the syringe.

Many designs of safety syringes have a noticeable gap between the zero graduation line and the position of the gasket which may cause immediate pre engagement before or during injection. There are a number of Disadvantages in using a fixed needle safety syringe with gap:

- The clinician cannot draw the medication from the zero graduation line.
- Difficult to aspirating precise medication, especially small volumes of medication.
- More difficult to remove the air bubbles from the syringe.
- Can not inject all the entire dose of medication without activating the safety mechanism.



**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316 Fax: +86 791-86496032  
E\_mail: info@kantecare.com MSN: kantecare@hotmail.com  
SKYPE: kantecare

[www.kantecare.com](http://www.kantecare.com)



**KANTECARE**  
we care, wherever you are



**Allergy Safety Syringe**    **Item No.:** KTC28-FNSS01

### Product Details

Kantecare Allergy safety syringe is a low-dose safety syringe, making it the perfect choice for the administration of allergy injections. Facilities using Kantecare Allergy Safety Syringe for allergy testing administration will be able to take advantage of the proven benefits of this safer alternative to traditional allergy syringe products.

Kantecare Allergy Tray is available in a tray of 25 1 ml or 0.5ml fixed needle Allergy safety syringes, making traditional repeated injection testing convenient.

### Special Product Features

- Special designed cap allow easy recap by using a one hand scoop technique to recap the needle (only if your facility allow recapping a clean needle for medical transport).
- Minimum waste volume to ensure accurate dosing.
- Allow low-angle intradermal injections and regular subcutaneous injections.
- Packed in a tray to help reduce waste for those facilities focusing on environmental friendly.
- Initiatives and make the traditional repeated allergy testing convenient.
- The safety feature engages with a single-handed technique.
- The clinician's hands and fingers remain behind the exposed sharp at all times.
- The clinician can easily tell whether the safety feature is activated.
- The needle is securely locked inside the safety device and remains protective through disposal.

**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: info@kantecare.com    MSN: kantecare@hotmail.com

SKYPE: kantecare

[www.kantecare.com](http://www.kantecare.com)







**KANTECARE**  
we care, wherever you are



**Insulin Safety Syringe**

**Item No.:** KTC28-FNSS02

### Product Details

Kantecare Fixed needle low dose *insulin safety syringes* is highly accurate and have little waste space volume, which leads to less medication waste. Upon completion of the insulin injection, the user easily renders the syringe safe by pulling back the plunger until it stops, at which point the needle is safely enclosed in the barrel. Then, snapping off the plunger to permanently destruct the plunger. Our *Insulin safety syringe* is available with 30U, 50U and 100U.

### Important Product Features

- Minimum waste volume to ensure accurate dosing.
- The safety feature engages with a single-handed technique.
- The clinician's hands and fingers remain behind the exposed sharp at all times.
- The clinician can easily tell whether the safety feature is activated.
- The needle is securely locked inside the safety device and remains protective through disposal.
- Orange needle cap clinically indicates a U-100 syringe.

**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: info@kantecare.com MSN: kantecare@hotmail.com

SKYPE: kantecare



[www.kantecare.com](http://www.kantecare.com)



**KANTECARE**  
we care, wherever you are



### Influenza Vaccine Safety Syringe

Item No.: KTC28-FNSS03

#### Product Details

Kantecare *Influenza Vaccine Safety syringe* with fixed needle designed specifically for institutions and companies that administer flu shots, such as long-term care facilities, health departments, and physician offices. Kantecare Influenza Vaccine Safety syringe with needle is available with both Trays which contain 25 1ml safety syringes and individual packing.

#### Important Product Features

- Minimum waste volume to ensure accurate dosing and reduce the waste of medication.
- Administer flu shots quickly and safely.
- The safety feature engages with a single-handed technique.
- The clinician's hands and fingers remain behind the exposed sharp at all times.
- The clinician can easily tell whether the safety feature is activated.
- The needle is securely locked inside the safety device and remains protective through disposal.

**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: info@kantecare.com MSN: kantecare@hotmail.com

SKYPE: kantecare



[www.kantecare.com](http://www.kantecare.com)



**KANTECARE**  
we care, wherever you are



**Tuberculin Safety Syringe**

**Item No.:** KTC28-FNSS04

#### **Product Details**

Kantecare fixed needle *tuberculin safety syringe* is the right device for tuberculin injection. Our fixed needle tuberculin safety syringe is available with 26 and 27g needles for Tuberculin injection.

#### **Special Product Features**

- Minimum waste volume is to ensure accurate dosing and reduce the waste of medication.
- The safety feature engages with a single-handed technique.
- The clinician's hands and fingers remain behind the exposed sharp at all times.
- The clinician can easily tell whether the safety feature is activated.
- The needle is securely locked inside the safety device and remains protective through disposal.

**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: info@kantecare.com MSN: kantecare@hotmail.com

SKYPE: kantecare



**www.kantecare.com**



**KANTECARE**  
we care, wherever you are



**Auto-Disable Syringe**

**Item No.:** KTC28-ADS

### Product Details

Our Auto-Disable syringe can prevent syringe reuse by locking and breaking the plunger after injection. The special safety feature makes it the right device for all kinds of injection, especial for large immunization program.

Available size: 0.5ml, 1ml, 2ml, 3ml, 5ml, 10ml

### Important Feature

- Permanently disabled.
- No add-on pieces.
- Looks and functions like a traditional syringe, No change in clinical technique.
- Hands and fingers stay behind the needle at all time.
- Accurate dose with minimal dead space.
- Bold, easy-to-read gradations.



**Kantech International (Jiangxi) Co., Ltd**

Tel: +86 791-86166316

Fax: +86 791-86496032

E\_mail: info@kantecare.com MSN: kantecare@hotmail.com

SKYPE: kantecare

[www.kantecare.com](http://www.kantecare.com)